



Rev. 6.4.2005

Code / Version: D007 / EN01

Clarity - Clarity Lite Comparison

Feature	Clarity	Clarity Lite
Support of multiple detector instruments	Up to 4 detectors / instrument	
Number of instruments connected simultaneously	1-4	1
Results summary table	yes	yes
Full history record including calibration data saved in a chromatogram	yes	yes
Time programmable integration parameters Peak Width and Threshold	yes	yes
Supporting tools for 21CFR Part 11:	yes	no
- Electronic signature	yes	no
- Audit trail	yes	no
- Advanced User Accounts	yes	no
Direct control (Autosamplers, GCs, HPLC pumps – <i>according to the List of Controlled Instruments</i>)	yes	no
System Suitability Module - SST	yes	no
Column Performance Table	yes	no
User Columns in Result table – Special calculation defined by the user	yes	no
Command Directories for selecting the working folder for data saving	yes	no
Archive/Restore functions	yes	no
Batch processing	yes	no
Commands Model method and Copy From for extended work with method files	yes	no
Compatibility with INT5 -ISA A/D converter	yes	no
Import/Export (AIA,TXT, dBase)	yes	yes
Postrun	yes	yes
Sequence	yes	yes
ISTD and ESTD calculation methods	yes	yes
Reference peaks usage	yes	yes
Automatical software update from Internet	yes	yes
Automatic subtraction of solvent defined in the method.	yes	no
Display background chromatogram in the Data Acquisition dialog	yes	no
Installation qualification (IQ) protocol	yes	no